



myActive Alert Mode Instruction Sheet

0.2
Rev

The Mytrex myActive Alert PHB has several modes to accommodate the needs of users. The steps required to enable each mode are shown below. Please refer to the myActive Alert User Manual for a description of each mode.

Supervised Fall Detection Mode

1. Press the myActive Alert button **5** times.
2. The LED on the myActive Alert will quickly flash green and red to verify that a mode was set.
3. **To verify that the mode was changed correctly:**
Press the myActive Alert button once. During the button press and 2 seconds thereafter the LED will show a fast, green blink.

Unsupervised Fall Detection Mode

1. This mode is set by pressing the myActive Alert button **7** times.
2. The LED on the myActive Alert will quickly flash green and red to verify that a mode was set.
3. **To verify that the mode was changed correctly:**
Press the myActive Alert button once. During the button press and 2 seconds thereafter the LED will show a steady and green.

Unsupervised RA3Tx Mode

1. This mode is set by pressing the myActive Alert button **10** times.
2. The LED on the myActive Alert will quickly flash green and red to verify that a mode was set.
3. **To verify that the mode was changed correctly:**
Press the myActive Alert button once. During the button press and 2 seconds thereafter the LED will show a slow, green blink.

Notes:

- *The amount of time allowed for every button press and time between button presses to change a mode must be less than ½ second or the mode change operation will cancel. It will also cancel a mode change operation if an invalid number of button presses is detected.*
- *The myActive Alert buttons are shipped in the “Unsupervised RA3Tx Mode”.*
- ***IMPORTANT! Ascertain that the mode was selected correctly by using the appropriate step 3 above. It is possible to program the wrong mode inadvertently if one or more button presses are not registered by the myActive Alert unit. (For example, thinking that 7 button presses have set the myActive Alert into Unsupervised Fall Detection Mode, when the unit only registered 5 button presses. Instead, this would set the mode to Supervised Fall Detection Mode, unbeknownst to the user, unless the step 3 verification procedure is performed.)***