

# MYTREX, INC.

SYNERGY BY DESIGN



## myActive Alert

The myActive Alert utilizes BioSensic's clinically academically validated ActivePERS technology, and proves to be the most accurate fall detection.

## Mytrex myActive Alert Features

- Compatible with Mytrex MXD, MXD3G and MXD-LTE
- Range of 600-1000 feet
- Transmits unique signals i.e. button press vs. fall-type event
- Auto-cancel feature (prevents a fall alert signal before to base unit) or (steps detected after will prevent signal from being sent to base unit)
- Monitored battery
- Can operate as a supervised Personal Help Button (PHB)
- Capable of sending unique alarm codes for fall-type events, button pushes, inactivity, supervisory alarms, and low battery alarms

**Note:** Due to the many types of falls that can occur, the user should always press the HELP button whenever help is needed and not rely solely upon the determination of the fall detection algorithms.

## Mytrex myActive Alert Features

- Dimensions: 2.165" x 1.468" x 0.603
- RF Transmission Frequency: 418 MHz
- Operating Temperature: -10°C to 60°C
- Battery Life: 1 year more (2 years in Basic non-supervised mode)
- Battery Voltage: 3.0V



877.801.6551

info@mytrexinc.com

