

myActive Alert:

The myActive Alert utilizes BioSensic's clinically and academically validated ActivePERS technology, and proves to be the most accurate fall detection

Pricing

○ 1-49.....	\$45.00
○ 50-99.....	\$42.00
○ 100-199.....	\$40.00
○ 200+.....	\$38.00

Specifications

- Dimension: 2.165" x 1.468" x 0.603"
- RF Transmission Frequency: 418 MHz
- Operating Temperature: -10°C to 60°C
- Battery Life: 1 year or more (2 years in Basic non-supervised mode)
- Battery Voltage: 3.0V

**Mytrex myActive Alert Features**

- Academically and clinically validated technology
- Can be used with the Mytrex MXD and MXD3G units
- Range of 600-1000 feet
- Send a fall detection and/or panic button signal
- Auto-cancel feature (prevents a fall alert signal before to base unit) or (Steps detected after will prevent signal from being sent to base unit)
- Monitored battery
- Can operate as a supervised Personal Help Button (PHB)
- Button Only Mode to preserve battery life
- Basic Mode
- Enhanced Mode
- Inactivity alert (Enhanced Mode)
- Capable of sending unique alarm codes for fall-type events, button pushes, inactivity, supervisory alarms, and low battery alarms
- Additional passive monitored features (Enhanced Mode)



877-801-6551

info@mytrexinc.com